



Health Professional Recommended
Medicare Code: E0628

Power Seat™

Keeps People Independent at Home

Uplift Power Seat™ is an electric-powered lift cushion that gently lifts and lowers users from an armchair or sofa.

TARGET CUSTOMER

- Customers who need help getting up due to arthritis, Parkinson's, MS, ALS or other joint and neuromuscular conditions.
- Customers who do not want to replace their favorite armchair or sofa.
- Customers who own a lift chair but need something to use in other rooms.
- People who are mostly confined indoors - whether at home, visiting, or in assisted living.

CUSTOMER BENEFITS

- Keeps people independent and at home longer.
- Can be stopped at any height needed.
- Conserves the users' energy for activities other than struggling to get up.

BEST SEATS TO USE

- A favorite armchair
- The sofa - beside the arm
- Dining room chair
- Kitchen chair

MEDICARE

- HCPCS Code: E0628, Medicare Part B
- Requires Certificate of Medical Necessity (CMN) from a physician.
- Conditions include arthritis, degenerative joint and neuromuscular disease.



Bringing quality to life.

UPLIFT

Power Seat™

FEATURES

EASY TO USE

The power lever acts like a switch. Simply click up to rise, push down to lower. The middle position will stop the lift at any time.



PATENTED ™



The flexible seat curves to a flat surface that lifts users up on a level plane. This patented design provides for a safe lift up, with no forward dumping motion.

POWER LEVERS



Power Seat comes with two detachable power levers that can be used on either the left or right side of the cushion.

WATERPROOF



The breathable fabric cover is waterproof and machine washable in cold water. Comes in black.

EASY CARE



The memory foam cushion has a water resistant coating that can be wiped down with mild soap and water.

BUILT-IN HANDLE NON-SKID BOTTOM



The base has a built-in handle and a non-slip material covering that prevents the cushion from shifting on fabric seating surfaces.

MEMORY FOAM CUSHION




Slow release memory foam provides additional comfort and prevention against pressure sores.



Model PS 1000

Lifts users up to 300lb

UPLIFT POWER SEAT SPECIFICATIONS

Model Number	PS1000
Medicare Code	E0628
Weight Capacity	300lb/136kg
Warranty	1 Year Limited
Weight	Product 12lb/5.4kg Shipping 14lb/6.3kg
Seat Dimensions	Depth 19in/48cm Width 16in/41cm
Seat Height	Open 10in/25cm Flat 2in/5cm
Electrical	Power Supply 120 Vac 60Hz 0.6A CSA & C/US Approved Cord Length 10ft/3m 
Patent Numbers	5,316,370 & 6,702,383

Uplift Technologies Inc.
#19-10 Morris Drive
Dartmouth, NS Canada B3B 1K8
www.up-lift.com