



PREFABRICATED SPINAL

STEALTH PRO™ LSO



SUPPORT METER



FIRM


For patients requiring firm support, the Stealth PRO combines the Stealth with a rigid panel to provide support from sacrococcygeal junction to T9 vertebra. The anterior opening design is ideal for patients requiring intracavitary pressure to reduce load on intervertebral discs and sagittal control. The PRO back panel attaches to the Stealth brace with our secure elastic suspender design, which ensures the same fit every time.

TYPICAL INDICATIONS

- Burst Fx
- Post-Op Spinal Lumbar Stabilization
- Post-Op Disk Surgery
- Chronic Back Pain
- Degenerative, Herniated or Bulging Disc
- Spondylolysis
- Spondylolisthesis
- Spinal Stenosis
- Facet Syndrome
- Tumors or Inflammation
- Spinal Instability
- Post-Op Laminectomy

SIZE	WAIST MEASUREMENT
SM	23" - 30"
MD	27" - 34"
LG	32" - 40"
XL	37" - 44"
2X	41" - 48"
3X	45" - 52"
4X	49" - 56"
MODEL #	DESCRIPTION
STPRO2W __	STEALTH PRO2 WHITE
STPRO4W __	STEALTH PRO4 WHITE
STPRO2B __	STEALTH PRO2 BLACK
STPRO4B __	STEALTH PRO4 BLACK

STANDARD COMPONENTS



CORSET



ANTERIOR POLYMER PANEL



PRO POSTERIOR POLYMER PANEL

RIGID PANELS

SEE P. 11 FOR STEALTH HEIGHT CHART